# BIODDT

#### Keywords: Protein Arrays, Manufacturing Process Development

#### Abstract

The development and optimization of manufacturable protein arrays As a manufacturer of nanoliter dispensing equipment requires researchers to address numerous technical challenges. The BioDot has had extensive experience with the developmen elements that contribute to the appropriate performance of arrays of manufacturing processes for the protein array industry. from a manufacturing perspective include the design of the array This poster highlights some of the main issues relating to itself, the type of dispensing system used, the environment in which appropriate process design for manufacturable protein arrays, dispensing is carried out and appropriate choices of substrates, and provides data demonstrating solutions to the issues reagents and buffers. Careful consideration must also be given to addressed, including data on the contribution of appropriate process design and process parameters such as cleaning of the system platform design, substrate choice, buffer choices, arrayer and degassing and filtering of reagents. As a result, optimizing the cleaning, system degassing, and environmental control, key process parameters to manufacture high quality protein arrays can be factors to the manufacturer's ability to robustly produce time consuming and costly.

protein microarrays.

#### 1. Array Design and Platform Scalability

#### Introduction

Commercially available nano- and pico-liter dispenser platforms issues. The BioJet Plus technology utilizes the typically are required to perform a wide variety of applications at a synchronization between the action of a high precision variety of scales, from laboratory to full scale production. Platform and pump, microsolenoid valve and robotic motion to create dispenser design must exhibit some degree of flexibility to overcome highly accurate and repeatable dispenses in the nanoliter these challenges.

In highly demanding arraying applications, dispensing platforms is adaptable to a variety of dispensing platforms and is typically should include a number of features critical to the process. completely scalable from a single channel on a small R&D These include aspirate and dispense capability, wash and dry stations platform to a multi-channel overhead gantry production for dispense tips carefully designed to avoid carry-over issues, humidity system without significant revalidation of the dispensing control, and the ability to customize deck configurations for different process. substrates.


The dispensing technology used is obviously critical to the success of into account, is critical to ensuring a smooth transition from the arraying process. BioDot's BioJet Plus technology is designed to research protein microarrays into a scaled up production line. provide a scalable solution to low nanoliter to microliter dispensing

range that can be placed on a substrate with extremely high positional accuracy at high speed. This dispensing technology


Careful process design, taking equipment design limitations



Flexibility of the design: BioDot AD5000 manufacturing printer, AD3200 research development printer and AD1500 research printer. BioDot arrayers are designed to accommodate different dispensing platforms for slides, microtiter plates, biochips, membranes, sheets of biosensors, etc...



Scalability of dispensing channels for faster throughput: BioDot's arrayer equipped with one dispensing channel can print twenty microarrays in 6 min 14 seconds, each array being 8 x 12 dots. It would take 3 min 7 seconds to do the same spotting experiment with two dispensing channels (14 nL spots, 1.02 mm center-to-center distance, 18 mm spacing between arrays).



Dot volume spot size:

500 µm

800 µm



1400 µm

#### **Experiment: Methods and Results**

Necessity to use different sizes of dispensing channel to obtain good protein microarrays over the 20 – 500 nL range of spot volumes. BioDot printer can accommodate different size of dispensing tips (the most common being 100 micron and 190 micron orifices). Larger spot size will suffer in their morphology if the size of the dispensing channel is not appropriate.

Arraying performed on BioDot AD5000 equipped with one dispensing channel and a humidity chamber. Eight different concentrations of pRagIgg were spotted under 50% humidity onto amine reactive OptArray slides. Print Buffer: 0.2M NaPO4 pH8.5, 0.005% CHAPS, 0.5% glycerol, and 0.01% BSA. Protein concentration doubled from row to row (bottom was 1.6 ng/mL top is 200 ng/mL). All slides were dried for two hours in a 50% humidity chamber. The slides were blocked for 1 hour in 50 mM NHEtOH and 1x PBS, then rinsed, then incubated with gIgg-Cy5 at 10 ng/mL in 1x PBS for 1 hour. The images were obtained on a Tecan LS400 scanner, at 633/670 nm, PMT 125, 10 microns.

Volume delivered (nL)	Average CV% over 960 drops	Stdev over 960 drops
20	7.10	3.8
30	7.52	4.0
100	5.00	2.6
500	1.41	0.6


Table 1: Coefficient of Variation (%CV) measured by fluorescence detection method on BioDot AD5000 equipped with 8 dispensing channel. Each %CV measures the variability in the volume delivered over 960 drops (except noted otherwise). Each dispensing channel delivers 120 drops.\* over 840 drops only

For higher density microarrays, a non-contact piezo-electric arrayer will deliver high throughput and high quality prints.



Robust technology for picoliter non contact dispensing: Arraying performed on BioDot SciFlex, 300 pL drop, 400 m center-to-center pitch. CVs are estimated to be around 1%. SciFlex robust technology is based on fluid pressure, piezo-electric impulse, and glass capillary.


clogging



	Tip Size 100 x 100µm	Tip Size 75 x 75μm	Tip Size 50 x 50μm
Average spot diameter	100-130 µm	75-110µm	50-80µm
% CV	7%	5%	4%
Total number of spots per one dip	350-400	400-500	500-600
Volumetric uptake	$\sim 0.1 \mu L$	$\sim 0.1 \mu L$	~0.1µL

sie 2 Statistical analysis of allays printed with SFF ins (Arrays are printed on one slide in test print mode)

### Hélène Citeau, Byron Doneen, Barbara McIntosh, Brendan O'Farrell BioDot, Inc. 17781 Sky Park Circle, Irvine, CA 92614 949.440.3685 www.biodot.com



By contrast, pin arrayers operate in contact mode and are thus more prone to

Contact-dispensing: Microarrays printed at PSTI using 75X75µm silicon pins spaced on 4.5mm centers with a spot spacing of 170µm (above). The above image shows all the spots (including prespotting) printed from a single uptake volume of 100nL. Arrays printed with a spot spacing of 145µm (below) using 75x75µm tips.

#### 2. Substrate Choice


#### Introduction

The substrate on which the protein microarrays will be printed is critical when surface than on another, causing a weaker signal, possible coalescence of evaluating a robotic arrayer for protein microarrays. Different glass slides carry the spots, and splattering. Identification of the best possible substrate for different surface chemistries. These distinct surface chemistries are responsible for a specific protein arraying process is critical to success. different interactions with the protein microarray: the spots may spread more on one

> inosilane gaps slide #1 Aminosilane gaps slide #3 0 0 0 -

#### Experiment: Methods and Results


**Illustration of different spot sizes and qualities in a protein microarray spotted** microarrays containing 9 dots of 20 nL each of protein solutions (125) on two different types of slides: epoxysilane coated glass slides and aminosilane ug/mL of goat anti-mouse in 0.05 M carbonate coupling buffer). Upon gaps glass slides. Arraying performed on BioDot AD3200 unit equipped with one developing, the slides were scanned on an Axon 4200AL scanner using dispensing channel and an in-line degasser. Each slide was spotted with fourteen laser settings of 40 and PMT of 230.




#### Fast Slide

Spot size differences between antibody arrays on S&S Fast Slides and on Teflon was used to change the refractive index of the membrane and to enhance poly-L-lysine coated slides. Arraying performed on BioDot AD5000 arrayer signal prior to scanning. Dot size was estimated to be between 400 and equipped with one dispensing channel. Each slide was spotted with 6x6 arrays of 12 550 microns on Fast Slides, and between 700 and 900 microns on Teflon nL dots of mouse HCG at 6 different concentrations in PBS buffer. Upon printing, coated poly-L-lysine slides. The protein concentration was 0, 0.12, 0.6, the slides were dried, blocked using 1% FCS in PBS, washed, incubated with 1 1.2, 6 and 12 ng from left to right on the array. ug/mL of goat anti-mouse IgG conjugated with AlexaFluor488, washed, dried and scanned on an IsoCyte<sup>™</sup> laser (488 nm). For optimum picture quality, mineral oil

3. Arraying Buffer Choice




#### Teflon poly-L-lysine coated slides



user must keep in mind that the manufacturer usually optimizes the buffer for one about 28% of the total diameter of the spot (see table 1). For all these type of glass substrate. Protein array wash buffer and protein array blocking buffer reasons, the user would greatly benefit from developing his own are also critical to the process. The latter is usually developed for one specific optimized buffer for his specific protein microarray. BioDot offers simple detection method. Protein array buffers are typically designed to enhance protein programs to achieve this iterative process efficiently using robotic stability and signal intensity. However, the importance of arraying buffers on the nanoliter dispensing equipment. spot roundness and size for protein microarrays is high: depending on the arraying

Several arraying buffers are commercially available for protein microarrays, but the buffer, in our hands spot diameter has been seen to vary by up to 200 um,



(5) PBS + 0.1% Tween-20

6) PBS + 0.1% SDS

Experiment: Methods and Results. Evaluation of seven spotting solutions. The spot Blocking Buffer, incubated with a biotinlated antibody derived from roundness and size for each arraying buffer was evaluated at 10 nL drop volumes. All goat, washed, incubated with streptavidin-Cy5 conjugate, washed and solutions contained the same concentration of antibodies (0.125 mg/mL of donkey dried. The slides were then scanned with P.E. Scan Array4000 imager. anti-goat antibody) but had different salts, detergents, sugars and other buffers The buffer composed of PBS, 10% glycerol, and 0.1% Triton X - 100 additives. The arrays were generated on BioDot AD5000 with one BioJetPlus vields spotted with the best combination of roundness, spot size dispensing channel equipped with a helium degasser. The dispensing height was consistency, and morphology. This buffer had the highest spot intensity, 0.14 mm from the surface of the slide. Upon printing, the slides were dried in the best aspect and heterogeneity measure, and was near the top in both desiccating chambers, shipped, and then blocked with 1 x S&S Protein Array

roundness and margination. However, it also yields the largest spots.

## **Design and development of manufacturable protein arrays**


Buffer #	Specific intensity	<b>CV intensity</b>	Ave Diameter (um)	STDEV	CV
(1)	40065	8	777	47	6.1%
(2)	35794	7	555	54	9.7%
(3)	26003	9	604	43	7.1%
(4)	26460	4	644	69	10.8%
(5)	23608	5	648	34	5.2%
(6)	32597	10	716	35	4.8%
(7)	37659	10	740	35	4.8%

Table 1: average spot diameter of protein in various print buffers


#### 4. Arrayer Cleaning -

#### Introduction

Mechanical cleaning of proteins from dispensing tips can be achieved using a costs and rate can be compromised in both contact and non-contact simple aspirate-and-dispense wash strategy. However, when defining a cleaning arraying when concentrated viscous reagents containing protein, strategy, the customer should consider the particular characteristics of their nucleic acid, binders and other polymers are used. In the contact proteins and choose cleaning reagents appropriately. Certain proteins or protein approach, production time allotted to pin cleaning can exceed the diluents may be incompatible with water and prefer organic based cleaning arraying period. Non-contact printing generally suffers fewer problems solutions. Customers should validate their own cleaning procedures for this with reagent clotting because liquid is exposed to the evaporative effects reason. In high throughput non-contact arraying, dispensing micro-droplets (nL of air for much briefer intervals. However, when highly concentrated to uL) of concentrated aqueous solutions of protein can cause blockage of solutions are used and hardware washing is infrequent, tip blockage ceramic tips with subsequent damage. Accuracy, precision and also production with subsequent over-pressurization induces print failures.



#### BioDot has developed and validated the efficiency of its own mechanical in-process wash condition, even after only a single wash cycle. IgG was used as a cleaning procedures to prevent carryover of protein between aspirated samples. In representative "sticky" protein. The alkaline phosphatase method of the array shown above, a full circle represents a 20 nL drop of antibody IG-AP, an development on Zeta Grip chips was chosen because of previously empty circle represent a 20 nL drop of buffer. Only spots in lines containing the demonstrated sensitivity and ease of use. The results indicate that the original IgG-AP reagent were visible upon development of the slides. No spots BioJet Plus ceramic tips do not retain substantial amounts of protein and developed in any of the lines containing wash solution after washing of the tips. simple water rinsing of the tips is sufficient to remove residual IgG after



Efficient washing steps: Elimination of protein by various washing agents when compared to water. Water by itself is capable of washing protein from the system. Study performed on BioDot AD5000.

Dependency of spot size on substrate. IgG in PBS on  $\zeta$ -Grip<sup>TM</sup>, (Miragene) and on Fast Slides (Schleicher and Schuell). Smaller spots on  $\zeta$ -Grip<sup>TM</sup> had a higher intensity, which allowed them to be read by a simple colorimetric method, instead of a





(7) PBS + 0.1% Triton X-100

#### 5. System Degassing

Experiment: Methods and Results

#### Introduction

Non-contact dispensing technology exhibits high level of reliability for continuous -meaning in the absence of any air bubbles. small volume dispensing under the proper experimental conditions. The De-aeration of the backing fluid is thus critical to maintain the accuracy and the precision of the drop volume are ensured by the steady-state pressure (SSP) inside the dispensing channel. This pressure will only be steady if the column of fluid inside the dispensing channel is




Flow-through vacuum degassing chamber with a single amorphous perfluorinated copolymer (Teflon<sup>®</sup> AF) degassing membrane. This compact design comprises continuously vented mini vacuum pump with a unitary PTFE diaphragm. Degassing efficiency is equal to fifty (50) times the degassing efficiency of other PTFE designs. Several vacuum degassing modules may be mounted in parallel on one arrayer. Eliminate any set-up time required for degassing prior to any dispensing experiment.

fluid path of non-contact dispensers free of any air bubble. The presence of any air bubble would lead to inaccurate dispensed volumes

#### Experiment: Methods and Results

Efficiency of vacuum chamber degassing




Efficiency of flow-through vacuum degassing chamber to lower the level of dissolved gas inside a fluid. Experiment performed with BioDot syringe pump. The fill flow rate is 30 uL /sec. The prime flow rate is 60 uL/sec. The vacuum efficiency is constant overtime and guarantees the degassing efficiency.

#### 6. Environmental Control

#### Introduction

During and after arraying, it is crucial to maintain a controlled deactivate the proteins and thus produce false negatives. The environment. Clean rooms are usually the best environment to obtain high addition of glycerol inside a protein solution may be useful to quality protein microarrays. The presence of dust during the arraying avoid the "doughnut effect" but it creates several challenges. Most process may cause splattering of the spots due to particulates on the commercially available arrayers have difficulties spotting solutions substrate surface or clogging of the dispensing nozzle. Once a protein array with high glycerol content. BioDot is capable of arraying high has been spotted successfully, the drying process deserves particular quality protein microarrays containing up to 10 % glycerol using attention. If this process is not thoroughly thought through, the protein BioJet Plus technology. However, the drying time is considerably array may turn out useless, as the formation of doughnuts instead of increased by the addition of glycerol. Unless the arraying buffer homogeneous spots will affect the interpretation of any quantization possesses glycerol, BioDot recommends leaving the slides spotted software. In most cases, it is necessary to maintain a relatively high with a protein array inside a humid environment for at least 1 humidity atmosphere (up to 60-70%) within the print enclosure to hour. This allows the protein array to dry slowly in a controlled prevent the spots from drying out from their external side. The drying out manner, which usually works well to avoid the "doughnut effect". of spots is deleterious in two ways: it creates a dry protein ridge called a "doughnut effect", which affects signal intensity measurements, and it may



#### **Experiment: Methods and Results**

Arraying performed on BioDot AD3200 equipped with one dispensing inside the humidity chamber for 1 hour after the completion of channel and flow-through vacuum degasser. The addition of 5-9% glycerol the array. Dot size was estimated to be around 800 microns. to their solutions helped avoiding drying out the 20 nL protein spot but Spots exhibited no "doughnut" effect in this experiment. the key to this issue resided mainly in the controlled drying environment. For this experiment, the humidity was set at 35% and the slides were left

			齫	· (III)
SI		ng		
(B))				
	line			

#### Incomplete drying process

#### Experiment: Methods and Results

Drying time is critical to the final quality of arrays spotted with glycerol. buffer. The size of each dot was estimated to be around 550 Arraying performed on BioDot AD5000 equipped with one dispensing microns. The humidity was set at 70% during the arraying to channel. One Nunc Maxisorp plate was arrayed with 10 nL dots of an avoid static effects. Imaging was performed by HRP conjugated commercial mouse antibody in a PBS/glycerol print

۲		

Good drying process (37C for 1 hr)

chemiluminescence after adding luminal to the plate.

### Conclusion

The effectiveness of any protein arraying experiment is dependent on the Plus non-contact, microsolenoid-based dispensing system, and interaction of a variety of components of the system, including the the Scienion piezo-electric based system can deposit dots of dispenser, the protein, the substrate, the binding chemistries/buffers, and 0.3-5000 nL, with CVs below 10%. the environmental and handling conditions.

The superiority of the printing method offered by BioDot lies in its consistency of printing for a wide range of heterogeneous samples. While contact printing methods may deposit dot volumes as low as 1 nL, with diameters in the order of 75 microns, contact dispensers are often limited by the viscosity of the solution and by clogging issues. BioDot's BioJet

In considering arraying systems, the user must keep in mind a variety of factors, as discussed above. Key among those considerations is the scalability of the manufacturing process, which is a reflection of the design and manufacturability of the product. Working with equipment manufacturers from an early stage of the product design is crucial to ensuring that process requirements and equipment capabilities match..